

anti- B4GALT6 antibody

Product Information

Catalog No.: FNab00776

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Required for the biosynthesis of glycosphingolipids.

Immunogen information

Immunogen: UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 6

Synonyms: None
Observed MW: 45-60kd
UniprotID: Q9UBX8

Application

Reactivity: Human, Mouse, Rat
Tested Application: ELISA, IHC, IF, WB

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human colon using FNab00776(B4GALT6 antibody) at dilution of 1:100



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00776(B4GALT6 Antibody) at dilution of 1:600