

# anti- BAG1L-speicfic antibody

## **Product Information**

Catalog No.:	FNab00787
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Rackground	

### Background

None

# **Immunogen information**

Immunogen:	BCL2-associated athanogene
Synonyms:	None
Observed MW:	None
Uniprot ID :	Q99933

## Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution	:IF: 1:20-1:200
Image:	

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)