

anti- BAIAP2L1 antibody

Product Information

Catalog No.:	FNab00794
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as adapter protein. Involved in the formation of clusters of actin bundles. Plays a role in the reorganization of the actin cytoskeleton in response to bacterial infection.

Immunogen information

Immunogen:	BAI1-associated protein 2-like 1
Synonyms:	BAIAP2L1, IRTKS
Observed MW:	57 kDa
UniprotID :	Q9UHR4

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab00794(BAIAP2L1 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00794(BAIAP2L1 Antibody) at dilution of 1:600