

# anti- BAMBI antibody

### **Product Information**

Catalog No.:	FNab00797
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

#### Background

Negatively regulates TGF-beta signaling.

### Immunogen information

Immunogen:	BMP and activin membrane-bound inhibitor homolog(Xenopuslaevis)
Synonyms:	BAMBI, Non metastatic gene A protein
Observed MW:	55kd
UniprotID :	Q13145

## Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution	:WB: 1:200-1:1000; IP: 1:200-1:1000
Image:	

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





IP Result of anti-BAMBI (IP:FNab00797, 3ug; Detection:FNab00797 1:300) with mouse spleen tissue lysate 4000ug.

rat kidney were subjected to SDS PAGE followed by western blot with FNab00797(BAMBI antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com