

anti- Barkor-Specific antibody

Product Information

Catalog No.:	FNab00806
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	KIAA0831
Synonyms:	ATG14, Atg14L, Barkor, KIAA0831
Observed MW:	55-60 kDa
Uniprot ID :	Q6ZNE5

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IF		
Recommended dilution: WB: 1:500-1:1000; IF: 1:10-1:100			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

1





2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com