

anti- BBOX1 antibody

Product Information

Catalog No.:	FNab00812
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the formation of L-carnitine from gamma-butyrobetaine.

Immunogen information

Immunogen:	butyrobetaine(gamma), 2-oxoglutarate dioxygenase(gamma- butyrobetaine hydroxylase) 1
Synonyms:	BBH, BBOX, BBOX1, G BBH, gamma BBH
Observed MW:	45 kDa
Uniprot ID :	O75936

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





human brain tissue were subjected to SDS PAGE followed by western blot with FNab00812(BBOX1 antibody) at dilution of 1:800

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com

2