

anti- BBS10 antibody

Product Information

Catalog No.: FNab00814

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Probable molecular chaperone. Assists the folding of proteins upon ATP hydrolysis. As part of the BBS/CCT complex may play a role in the assembly of BBSome, a complex involved in ciliogenesis regulating transports vesicles to the cilia. Involved in adipogenic differentiation.

Immunogen information

Immunogen: Bardet-Biedl syndrome 10

Synonyms: Bardet Biedl syndrome 10, BBS10, C12orf58

Observed MW: 81kd

Uniprot ID: Q8TAM1

Application

Reactivity: Human, Zebrafish

Tested Application: ELISA, WB, IF

Recommended dilution: WB: 1:500-1:1000; IF: 1:20-1:200

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1





HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00814(BBS10 antibody) at dilution of 1:800