

anti- BCRP/ABCG2 antibody

Product Information

Catalog No.:	FNab00855
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

Immunogen information

Immunogen:	ATP-binding cassette, sub-family G(WHITE), member 2
Synonyms:	None
Observed MW:	60-70 kDa
UniprotID :	Q9UNQ0

Application


Reactivity:	Human, Mouse		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00855(BCRP,ABCG2 Antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00855(BCRP,ABCG2 antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>