

anti- BCS1L antibody

Product Information

Catalog No.: FNab00857

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 9G10 IsoType: IgG2a

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Chaperone necessary for the assembly of mitochondrial respiratory chain complex III. Plays an important role in the maintenance of mitochondrial tubular networks, respiratory chain assembly and formation of the LETM1 complex.

Immunogen information

Immunogen: BCS1-like(yeast)

Synonyms: BCS1
Observed MW: 47 kDa
UniprotID: Q9Y276

Application

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1





Immunohistochemistry of paraffin-embedded human gliomas using FNab00857(BCS1L antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab00857(BCS1L antibody) at dilution of 1:1000