

anti- BEST1 antibody

Product Information

Catalog No.:	FNab00864
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Forms calcium-sensitive chloride channels. Highly permeable to bicarbonate.

Immunogen information

Immunogen:	bestrophin 1
Synonyms:	VMD2
Observed MW:	68 kDa
Uniprot ID :	O76090

Application


Reactivity:	Human		
Tested Application:	ELISA, WB		
Recommended dilution: WB: 1:500-1:2000			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00864(BEST1 antibody) at dilution of 1:1000



2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com