

anti- BHMT2 antibody

Product Information

Catalog No.:	FNab00891
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the regulation of homocysteine metabolism. Converts homocysteine to methionine using S-methylmethionine(SMM) as a methyl donor.

Immunogen information

Immunogen:	betaine-homocysteine methyltransferase 2
Synonyms:	BHMT2
Observed MW:	40-45 kDa
Uniprot ID :	Q9H2M3

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB		
Recommended dilution: WB: 1:500-1:2000			
Image:			


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00890(BHMT2 antibody) at dilution of 1:2000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)