

anti- BIRC8 antibody

Product Information

Catalog No.:	FNab00901
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Protects against apoptosis mediated by BAX.

Immunogen information

Immunogen:	baculoviral IAP repeat-containing 8
Synonyms:	BIRC8, hILP2, IAP like protein 2, ILP 2, ILP2
Observed MW:	27 kDa
Uniprot ID :	Q96P09

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00901(BIRC8 antibody) at dilution of 1:1000