

anti- BL21(DE3)-lysate antibody

Product Information

Catalog No.: FNab00903

Size: 100μg Form: liquid

Purification: Protein A purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: BL21(DE3)-lysate

Synonyms: BL21(DE3) lysate

Observed MW: 1-300 kDa

Application

Reactivity: Bacterial

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>