

# anti- BMP7 antibody

# **Product Information**

Catalog No.:	FNab00920
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

# Background

Induces cartilage and bone formation. May be the osteoinductive factor responsible for the phenomenon of epithelial osteogenesis. Plays a role in calcium regulation and bone homeostasis.

#### Immunogen information

Immunogen:	bone morphogenetic protein 7
Synonyms:	BMP 7, BMP7, bone morphogenetic protein 7, Eptotermin alfa, OP 1, OP1, Osteogenic protein 1
Observed MW:	49 kDa
Uniprot ID :	P18075

# Application


Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, IHC, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan,<br/>Hubei, China(430206)Tel :( 0086)027-87384275Fax: (0086)027-87800889www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human renal cell carcinoma tissue slide using FNab00920(BMP7 Antibody) at dilution of 1:50

Immunofluorescent analysis of HeLa cells using FNab00920 (BMP7 antibody) at dilution of 1:50 and Alexa Fluor 488-congugated AffiniPure Goat Anti-Rabbit IgG(H+L)

IP Result of anti-BMP7 (IP: FNab00920, 4ug; Detection: FNab00920 1:800) with HEK-293 cells lysate 2800ug.

HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00920(BMP7 Antibody) at dilution of 1:1000

### 2

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fa

Fax: (0086)027-87800889

www.fn-test.com