

anti- BTNL2 antibody

Product Information

| Catalog No.: | FNab00985 |
|---------------|---|
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | \geq 95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Negative regulator of T-cell proliferation.

Immunogen information

| Immunogen: | butyrophilin-like 2 |
|--------------|---|
| Synonyms: | BTL II, BTNL2, Butyrophilin like protein 2, HSBLMHC1, SS2 |
| Observed MW: | 65 kDa |
| Uniprot ID : | Q9UIR0 |

Application

| Reactivity: | Human, Mouse |
|-----------------------|------------------|
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:500-1:1000 |
| Image: | |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00985(BTNL2 Antibody) at dilution of 1:600

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)