

# anti- BubR1 antibody

#### **Product Information**

Catalog No.: FNab00989

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

## **Immunogen information**

Immunogen: budding uninhibited by benzimidazoles 1 homolog beta(yeast)

Synonyms: None

Observed MW: 120-130 kDa

UniprotID: O60566

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab00989(BUB1B Antibody) at dilution of 1:200



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00989(BUB1B antibody) at dilution of 1:500