

# anti- C10orf67 antibody

### **Product Information**

Catalog No.:	FNab01000
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

#### Background

None

## **Immunogen information**

Immunogen:	chromosome 10 open reading frame 67
Synonyms:	LINC01552
Observed MW:	20 kDa
UniprotID :	Q8IYJ2

## Application


Reactivity:	Human			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:200-1:1000				
Image:				

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01000( C10orf67 Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>