

anti- C20orf149 antibody

Product Information

Catalog No.:	FNab01080
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	chromosome 20 open reading frame 149
Synonyms:	C20orf149, EXPDF, PPDF, PPDF,EXPDF
Observed MW:	12 kd
Uniprot ID :	Q9H3Y8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon using FNab01080(C20orf149 antibody) at dilution of 1:100



Immunofluorescent analysis of HeLa cells, using FNab01080(C20orf149 antibody) at 1:25 dilution. Blue pseudocolor = DAPI (fluorescent DNA dye).



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab01080(C20orf149 antibody) at dilution of 1:300