

anti- CACNG1 antibody

Product Information

Catalog No.:	FNab01177
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This protein is a subunit of the dihydropyridine(DHP) sensitive calcium channel. Plays a role in excitation-contraction coupling. The skeletal muscle DHP-sensitive Ca(2+) channel may function only as a multiple subunit complex.

Immunogen information

Immunogen:	calcium channel, voltage-dependent, gamma subunit 1
Synonyms:	CACNLG
Observed MW:	None
Uniprot ID :	Q06432

Application

Reactivity:	Human
Tested Application:	ELISA, IF, IHC
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01177(anti- CACNG1 antibody) at dilution of 1:50