

anti- CADM4 antibody

Product Information

Catalog No.:	FNab01192
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the cell-cell adhesion. Has calcium-and magnesium-independent cell-cell adhesion activity. May have tumor-suppressor activity.

Immunogen information

Immunogen:	cell adhesion molecule 4
Synonyms:	IGSF4C, NECL4, TSLL2
Observed MW:	65-70 kDa
UniprotID :	Q8NFZ8

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab01192(CADM4 antibody) at dilution of 1:1000