

anti- calreticulin antibody

Product Information

Catalog No.:	FNab01224
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	calreticulin
Synonyms:	CALR, Calregulin, calreticulin, cC1qR, CRP55, CRT, CRTC, ERp60, grp60, HACBP, RO, SSA
Observed MW:	60 kDa
Uniprot ID :	P27797

Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human thyroid tissue slide using FNab01224 (calreticulin Antibody) at dilution of 1:200



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab01224 (calreticulin Antibody) at dilution of 1:600