

anti- Casein Kinase 1 delta antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab01279 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|------------------------|
| Immunogen: | casein kinase 1, delta |
| Synonyms: | None |
| Observed MW: | 47kd |
| Uniprot ID : | P48730 |

Application


| | |
|-----------------------|----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IF, IP |
| Recommended dilution: | WB: 1:500-1:3000; IF: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-Casein Kinase 1 delta
(IP:FNab01279, 3ug; Detection:FNab01279
1:1000) with HeLa cells lysate 2000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab01279(Casein Kinase 1
delta antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com