

anti- CBLN2 antibody

Product Information

Catalog No.:	FNab01321
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play role in synaptogenesis induction.

Immunogen information

Immunogen:	cerebellin 2 precursor
Synonyms:	CBLN2, Cerebellin 2, cerebellin 2 precursor
Observed MW:	24 kDa
Uniprot ID :	Q8IUUK8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01321(CBLN2 antibody) at dilution of 1:800