

## anti- CCR5 antibody

### Product Information

Catalog No.:	FNab01391
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Receptor for a number of inflammatory CC-chemokines including MIP-1-alpha, MIP-1-beta and RANTES and subsequently transduces a signal by increasing the intracellular calcium ion level. May play a role in the control of granulocytic lineage proliferation or differentiation. Acts as a coreceptor(CD4 being the primary receptor) for HIV-1 R5 isolates. (Microbial infection) Acts as a receptor for human immunodeficiency virus-1/HIV-1.

### Immunogen information

Immunogen:	chemokine(C-C motif) receptor 5
Synonyms:	CMKBR5
Observed MW:	45 kDa
UniprotID :	P51681

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, FC
Recommended dilution:WB:	1:500-1:2000


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01391(CCR5 antibody) at dilution of 1:800

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)