

anti- CD98 antibody

Product Information

Catalog No.:	FNab01511
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	solute carrier family 3(activators of dibasic and neutral amino acid transport), member 2
Synonyms:	None
Observed MW:	85 kDa, 120 kDa
Uniprot ID :	P08195

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01511(SLC3A2 Antibody) at dilution of 1:500



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01511(SLC3A2 antibody) at dilution of 1:500