

## anti- Chromogranin A antibody

### Product Information

Catalog No.:	FNab01687
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6F9
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	chromogranin A(parathyroid secretory protein 1)
Synonyms:	CGA, CHGA, Chromogranin A, Pituitary secretory protein I, SP I, Vasostatin I Vasostatin 2
Observed MW:	50 kDa, 70 kDa
Uniprot ID :	P10645

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:1000-1:2000; IHC: 1:10000-1:40000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab01687(CHGA Antibody) at dilution of 1:10000



rat pancreas tissue were subjected to SDS PAGE followed by western blot with FNab01687(CHGA Antibody) at dilution of 1:2000