

## anti- Citrate synthase antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab01722   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

None

### Immunogen information

|              |                  |
|--------------|------------------|
| Immunogen:   | citrate synthase |
| Synonyms:    | None             |
| Observed MW: | 50 kDa           |
| UniprotID :  | O75390           |

### Application


|                       |   |
|-----------------------|---|
| Reactivity:           | Human, Mouse, Rat   |
| Tested Application:   | ELISA, WB, IHC, IF, IP  |
| Recommended dilution: | WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200 |
| Image:                |   |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01722(CS Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



IP Result of anti-Citrate synthase (IP:FNab01722, 4ug; Detection:FNab01722 1:1000) with mouse heart tissue lysate 3200ug.



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01722(CS antibody) at dilution of 1:2000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)