

anti- Claudin 18 antibody

Product Information

Catalog No.:	FNab01735
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3D2
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	claudin 18
Synonyms:	claudin 18, CLDN18, SFTA5, SFTPJ
Observed MW:	28 kDa
Uniprot ID :	P56856

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IP
Recommended dilution:	IHC: 1:200-1:1000; IP: 1:200-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab01735 (Claudin 18 antibody) at dilution of 1:600



IP Result of anti-Claudin 18 (IP: FNab01735, 3ug; Detection: FNab01735 1:500) with mouse stomach tissue lysate 3000ug.