

anti- CLOCK antibody

Product Information

Catalog No.:	FNab01769
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	clock homolog(Mouse)
Synonyms:	BHLHE8, KIAA0334
Observed MW:	110 kDa
Uniprot ID :	O15516

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP result of anti-CLOCK(IP:FNab01769, 4ug;
Detection:FNab01769 1:300) with HeLa cells
lysate 1600 ug



Various lysates were subjected to SDS PAGE
followed by western blot with
FNab01769(CLOCK antibody) at dilution of 1:500