
Cyclin D2 antibody

Product Information


Catalog No.:	FNab02125
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	cyclin D2
Synonyms:	CCND2, cyclin D2, G1/S-specific cyclin-D2
Observed MW:	36 kDa
Uniprot ID :	P30279

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded Skeletal muscle tissue using FNab02125(Cyclin D2 antibody) at dilution of 1:100



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab02125(Cyclin D2 antibody) at dilution of 1:600