

anti- CYP19A1 antibody

Product Information

Catalog No.:	FNab02146
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the formation of aromatic C18 estrogens from C19 androgens.

Immunogen information

Immunogen:	cytochrome P450, family 19, subfamily A, polypeptide 1
Synonyms:	ARO1, CYAR, CYP19
Observed MW:	49-55 kDa
UniprotID :	P11511

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab02146(CYP19A1 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02146(CYP19A1 antibody) at dilution of 1:1000