

## anti- Cytokeratin 15 antibody

### Product Information

Catalog No.:	FNab02199
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6E6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	keratin 15
Synonyms:	CK 15, CK15, Cytokeratin 15, K15, K1CO, keratin 15, KRT15, KRTB
Observed MW:	49 kDa
Uniprot ID :	P19012

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



A431 cells were subjected to SDS PAGE followed by western blot with FNab02199(Cytokeratin 15 antibody) at dilution of 1:1000