

anti- DAT antibody

Product Information

Catalog No.:	FNab02247
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Amine transporter. Terminates the action of dopamine by its high affinity sodium-dependent reuptake into presynaptic terminals.

Immunogen information

Immunogen:	solute carrier family 6(neurotransmitter transporter, dopamine), member 3
Synonyms:	DA transporter, DAT, DAT1, SLC6A3
Observed MW:	68 kDa
Uniprot ID :	Q01959

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab02247(SLC6A3 Antibody) at dilution of 1:50



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab02247(SLC6A3 Antibody) at dilution of 1:600