

anti- DCLK1 antibody

Product Information

Catalog No.:	FNab02266
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Probable kinase that may be involved in a calcium-signaling pathway controlling neuronal migration in the developing brain. May also participate in functions of the mature nervous system.

Immunogen information

Immunogen:	doublecortin-like kinase 1
Synonyms:	DCAMKL1, DCDC3A, KIAA0369
Observed MW:	48 kDa, 82 kDa
UniprotID :	O15075

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain using FNab02266(DCLK1 antibody) at dilution of 1:100



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab02266(DCLK1 Antibody) at dilution of 1:1000