

anti- DCTN4 antibody

Product Information

Catalog No.:	FNab02274
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Could have a dual role in dynein targeting and in ACTR1A/Arp1 subunit of dynactin pointed-end capping. Could be involved in ACTR1A pointed-end binding and in additional roles in linking dynein and dynactin to the cortical cytoskeleton.

Immunogen information

Immunogen:	dynactin 4(p62)
Synonyms:	None
Observed MW:	62 kDa
UniprotID :	Q9UJW0

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-DCTN4 (IP:FNab02274, 4ug;
Detection:FNab02274 1:700) with Jurkat cells
lysate 2800ug.



human brain tissue were subjected to SDS PAGE
followed by western blot with FNab02274(DCTN4
antibody) at dilution of 1:800