

anti- DDX3Y antibody

Product Information

Catalog No.:	FNab02307
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Probable ATP-dependent RNA helicase. May play a role in spermatogenesis.

Immunogen information

Immunogen:	DEAD(Asp-Glu-Ala-Asp) box polypeptide 3, Y-linked
Synonyms:	DBY
Observed MW:	73 kDa
Uniprot ID :	O15523

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis using FNab02307(DDX3Y antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02307(DDX3Y Antibody) at dilution of 1:300