

# anti- Desmocollin 2 antibody

#### **Product Information**

Catalog No.: FNab02344

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 7D7
IsoType: IgG2a

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen: desmocollin 2

Synonyms: ARVD11, Cadherin family member 2, CDHF2, desmocollin 2,

Desmocollin 3, Desmosomal glycoprotein II, Desmosomal glycoprotein

III, DG2, DGII/III, DSC2, DSC3

Observed MW: 110 kDa Uniprot ID: Q02487

#### **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB

Recommended dilution: WB: 1:1000-1:4000

Image:


#### 1

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02344(Desmocollin 2 antibody) at dilution of 1:2000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>