

## anti- DHX38 antibody

### Product Information

Catalog No.:	FNab02381
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Probable ATP-binding RNA helicase involved in pre-mRNA splicing.

### Immunogen information

Immunogen:	DEAH(Asp-Glu-Ala-His) box polypeptide 38
Synonyms:	DDX38, KIAA0224, PRP16
Observed MW:	141 kDa
UniprotID :	Q92620

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

IP Result of anti-DHX38 (IP:FNab02381, 4ug;  
Detection:FNab02381 1:700) with HeLa cells  
lysate 1400ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab02381(DHX38  
antibody) at dilution of 1:1000

