

anti- DNAJB2 antibody

Product Information

Catalog No.:	FNab02447
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily B, member 2
Synonyms:	HSJ1, HSPF3
Observed MW:	36-40 kDa
UniprotID :	P25686

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-DNAJB2 (IP:FNab02447, 4ug;
Detection:FNab02447 1:1000) with mouse brain
tissue lysate 4000ug.



K-562 cells were subjected to SDS PAGE
followed by western blot with
FNab02447(DNAJB2 antibody) at dilution of
1:300