

anti- DNAJC21 antibody

Product Information

Catalog No.:	FNab02464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May act as a co-chaperone for HSP70.

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 21
Synonyms:	DNAJA5
Observed MW:	70 kDa
UniprotID :	Q5F1R6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab02464(DNAJC21 Antibody) at dilution of 1:50



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab02464(DNAJC21 antibody) at dilution of 1:2000