

anti- DNASE2B antibody

Product Information

Catalog No.:	FNab02478
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Hydrolyzes DNA under acidic conditions. Does not require divalent cations for activity. Participates in the degradation of nuclear DNA during lens cell differentiation.

Immunogen information

Immunogen:	deoxyribonuclease II beta
Synonyms:	DLAD
Observed MW:	42 kDa
UniprotID :	Q8WZ79

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab02478(DNASE2B Antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02478(DNASE2B Antibody) at dilution of 1:2000