

anti- DPCD antibody

Product Information

Catalog No.:	FNab02512
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in the formation or function of ciliated cells.

Immunogen information

Immunogen:	deleted in a Mouse model of primary ciliary dyskinesia
Synonyms:	None
Observed MW:	23-26 kDa
UniprotID :	Q9BVM2

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart slide using FNab02512(DPCD Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab02512(RP11-529I10.4 antibody) at dilution of 1:300