

anti- DPP7 antibody

Product Information

Catalog No.:	FNab02521
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays an important role in the degradation of some oligopeptides.

Immunogen information

Immunogen:	dipeptidyl-peptidase 7
Synonyms:	DPP2, QPP
Observed MW:	54kd
UniprotID :	Q9UHL4

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02521(DPP7 antibody) at dilution of 1:300