

anti-Emerin antibody

Product Information

Catalog No.:	FNab02753
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	emerin
Synonyms:	None
Observed MW:	35kd
UniprotID :	P50402

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP result of EMD antibody (FNab02753 for IP and Detection) with mouse skeletal muscle tissue lysate.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02753 (Emerin antibody) at dilution of 1:1000