

## anti- ERC2 antibody

### Product Information

Catalog No.:	FNab02831
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Thought to be involved in the organization of the cytomatrix at the nerve terminals active zone(CAZ) which regulates neurotransmitter release. Seems to act together with BSN.May recruit liprin-alpha proteins to the CAZ.

### Immunogen information

Immunogen:	ELKS/RAB6-interacting/CAST family member 2
Synonyms:	KIAA0378
Observed MW:	115-120kd
UniprotID :	O15083

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02831(ERC2 antibody) at dilution of 1:500