

## anti- E-selectin /CD62E antibody

### Product Information

Catalog No.:	FNab02863
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	selectin E
Synonyms:	CD62E, ELAM, ELAM 1, ELAM1, ESEL, E-selectin, LECAM2, SELE, selectin E
Observed MW:	67 kDa
Uniprot ID :	P16581

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:200-1:2000; IHC: 1:50 - 1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded Human kidney using FNab02863(E-selectin antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02863(E-selectin antibody) at dilution of 1:300