

anti- Evi2a antibody

Product Information

Catalog No.:	FNab02885
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May complex with itself or/and other proteins within the membrane, to function as part of a cell-surface receptor.

Immunogen information

Immunogen:	ecotropic viral integration site 2a
Synonyms:	Evi-2, Evi-2a, Evi2
Observed MW:	26 kDa
UniprotID :	P22794

Application


Reactivity:	Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02885(Evi2a antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com