

## anti- F5 antibody

### Product Information

Catalog No.:	FNab02923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Central regulator of hemostasis. It serves as a critical cofactor for the prothrombinase activity of factor X that results in the activation of prothrombin to thrombin.

### Immunogen information

Immunogen:	coagulation factor V (proaccelerin, labile factor)
Synonyms:	None
Observed MW:	145 kDa
UniprotID :	P12259

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02923(F5 antibody) at dilution of 1:1000