

anti- FABP3 antibody

Product Information

Catalog No.:	FNab02931
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

Immunogen information

Immunogen:	fatty acid binding protein 3, muscle and heart(mammary-derived growth inhibitor)
Synonyms:	FABP11, FABP3, Fatty acid binding protein 3, H FABP, HFABP, M FABP, MDGI, O FABP
Observed MW:	15 kDa
Uniprot ID :	P05413

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02931(FABP3 Antibody) at dilution of 1:1000