

anti- FAF1 antibody

Product Information

Catalog No.:	FNab02946
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9H10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Potentiates but cannot initiate FAS-induced apoptosis.

Immunogen information

Immunogen:	Fas(TNFRSF6) associated factor 1
Synonyms:	UBXD12, UBXN3A
Observed MW:	74kd
Uniprot ID :	Q9UNN5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02946(FAF1 antibody) at dilution of 1:1000

